



BYRNA CL HOLSTER

YOUR BYRNA CL HOLSTER COMES WITH:

- (1) Pair of IWB Belt Loop Attachment (for Inside the Waistband)
- (1) Pair of OWB Belt Loop Attachment (for Outside the Waistband)
- Adjustable Retention Mechanism

RETENTION ADJUSTMENT INSTRUCTIONS:

LOCATE RETENTION SCREWS:

The two retention screws are located near the trigger guard area on the holster.

ADJUST FOR DESIRED RETENTION:

- **Increase Retention:** Turn both screws clockwise with the Allen key to tighten the grip on your Byrna launcher.
- **Decrease Retention:** Turn both screws counter-clockwise with the Allen key to loosen the grip.

TEST THE DRAW:

After adjusting, insert and draw your launcher multiple times to ensure the holster holds it securely while allowing a smooth draw.

NOTE:

The Holster can be set to be Left or Right handed by reversing the belt loop attachments.



BYRNA CL HOLSTER

**Scan the QR code for more
Instruction and Videos.**

